



The ResQGARD ITD is used to treat low blood circulation, which results from hypotension. Hypotension has a significantly detrimental impact on morbidity and mortality. Up to 20% of patients admitted to the hospital from the emergency department (ED) have a hypotensive episode, and **patients with hypotension are ten times more likely to have sudden and unexpected death.**²² A large national trauma bank study found that for every 10 mmHg decrease in systolic blood pressure (BP), mortality increased 4.8%.²⁸

The ResQGARD is an impedance threshold device (ITD) that harnesses the body's natural reflexes to enhance circulation. Over 40 published clinical trials, animal studies and review articles, many in collaboration with NASA and/or the US Department of Defense, have shown that the ResQGARD works to rapidly and noninvasively increase blood pressure from a variety of etiologies (e.g. hypovolemia, orthostatic intolerance, heat shock, renal dialysis and blood donation). In 2008, the ResQGARD received the US Army's SBIR Achievement Award for the technology's application in the non-invasive treatment of hemorrhagic shock.

Human Clinical Trials

The ResQGARD has been evaluated in 20 human clinical trials including models of normotension (6,7,18,24), orthostatic intolerance (6,7,18,24,25,36,38) and hypotension:

- In an ED and emergency medical services (EMS) setting (20,34,36,38,40,41)
- During simulated hemorrhagic shock (23,27,30)
- During orthostatic challenge (e.g. supine to standing or tilt-table/squat-stand test) (8,11,12,15,17,25,29)
- During renal dialysis (14) and pregnancy (41)
- Following blood donation (14)

These studies have shown that the ResQGARD:

- Improves hemodynamics
 - Increases systolic and diastolic BP 4 - 30% (4,6,7,14,18,20,23,36,38,40,41)
 - Increases mean arterial pressure (MAP) 5 - 27% (4,15,18,23,34,36,38)
 - Increases stroke volume and cardiac output 10 - 21% (7,15)
 - Reduces the drops in pressure, cardiac output and/or stroke volume during orthostatic stress (11,12,15,17,25,27,30)
 - Increases cerebral blood flow 9 - 11% (18,27,29)
- Increases BP in actual or simulated hemorrhagic shock, but not to levels typically associated with "popping the clot" (13,14,20,23,27,30,32,36,37,38,40)
- Reduces orthostatic symptoms (15,25,29)
- Lowers intrathoracic pressure during inspiration (6,7,18,23,24)
- Breathing through therapeutic resistance is well tolerated (6,7,14,23,24,36,38,40)
- Does not compromise oxygen saturation (23,38) or cause increase in heart rate (36)

Animal Studies

The ResQGARD has been evaluated in 12 animal studies (1,3,4,9,10,19,21,26,31,34,35,39) including animal models of hypovolemia and hemorrhagic shock (1,4,9,19,21,35,39), pediatric models of hypovolemia (4[w/PEEP],26), and models of heat shock (31).

These studies have shown that the ResQGARD:

- Lowers intrathoracic pressures during inspiration (3,4,9,10,19,21,31)
- Improves hemodynamics
 - Increases systolic and diastolic BP (1,4,9,10,31,34,35,39)
 - Increases MAP (4,19,31,34,35,39)
 - Lowers intracranial pressure during inspiration (3,21)
 - Increases cerebral perfusion pressure (3,21)
 - Increases coronary perfusion pressure (4,31)
 - Lowers right atrial pressure (3,21)
 - Increases cardiac output/index (9,10,34)
 - Increases stroke volume index (10)
 - Increases end tidal carbon dioxide (1,9,31)
 - Does not compromise arterial oxygenation (19,31)
- Compared to fluid therapy, does not dilute clotting factors (39)
- Increases blood pressure in hemorrhagic shock, but not to levels typically associated with "popping the clot" (9,10,19,21,26,35,39)
- Increases work of breathing but is well tolerated (26)
- "Buys time" and extends the "window of opportunity" during shock until fluids can be administered (1,9,10,19,31,39)
- Increases short-term (90-min) (19,26) and 24-hour survival (26)

Bibliography

1. **Lurie KG** et al. Use of an inspiratory impedance threshold valve for the rapid treatment of hemorrhagic shock in spontaneously breathing pigs. *Crit Care Med* 2001;29(12)Suppl:A13.
2. **Lurie KG** et al. Augmentation of ventricular preload during treatment of cardiovascular collapse and cardiac arrest. *Crit Care Med* 2002;30(4 Suppl):S162-165.
3. **Lurie KG** et al. Cardio-cranial interactions: reduction of elevated intracranial pressure with an inspiratory ITD in spontaneously breathing pigs after resuscitation from cardiac arrest. *Crit Care Med* 2003;31(12)Suppl:A10.
4. **Samniah N** et al. Feasibility and effects of transcutaneous phrenic nerve stimulation combined with an inspiratory impedance threshold in a pig model of hemorrhagic shock. *Crit Care Med* 2003;31(4):1197-1202.
5. **Sigurdsson G** et al. Cardiorespiratory interactions and blood flow generation during cardiac arrest and other states of low blood flow. *Curr Opin Crit Care* 2003;9:183-188.
6. **Convertino VA** et al. Effects of inspiratory impedance on the carotid-cardiac baroreflex response in humans. *Clin Auton Res* 2004;14:240-8.
7. **Convertino VA** et al. Hemodynamics associated with breathing through an inspiratory ITD in human volunteers. *Crit Care Med* 2004;32(9):S381-366.
8. **Doerr DF** et al. Comparison of two methods of non-invasive determination of cardiac output during an orthostatic challenge. *Aviat Space Environ Med* 2004;75:B117.
9. **Lurie KG** et al. Treatment of hypotension in pigs with an inspiratory ITD: a feasibility study. *Crit Care Med* 2004;32(7):1555-1562.
10. **Marino BS** et al. Spontaneous breathing through an inspiratory ITD augments cardiac index and stroke volume index in a pediatric porcine model of hemorrhagic hypovolemia. *Crit Care Med* 2004;32(9 Suppl):S398-405.
11. **Melby DP** et al. An inspiratory impedance device improves hemodynamics during a squat-stand test: implications for treating orthostatic intolerance. *Heart Rhythm* 2004;1(18):S257.
12. **Melby DP** et al. A novel inspiratory ITD may diminish orthostatic intolerance. *Heart Rhythm* 2004;1(18):S227.
13. **Convertino VA** et al. Inspiratory resistance as a potential treatment for orthostatic intolerance and hemorrhagic shock. *Aviat Space Environ Med* 2005;76(4):319-325.
14. **Convertino VA** et al. Restoration of central blood volume: application of a simple concept and simple device to counteract cardiovascular instability in syncope and hemorrhage. *J Gravit Physio* 2005;12(1):P55-60.
15. **Convertino VA** et al. Effects of inspiratory impedance on hemodynamic responses to a squat-stand test in human volunteers: implications for treatment of orthostatic hypotension. *Eur J Appl Physiol* 2005;94(4):392-399.
16. **Aufderheide TP** et al. Vital organ blood flow with the ITD. *Crit Care Med* 2006;34(12 Suppl):S466-473.
17. **Convertino VA** et al. Inspiratory impedance effects on hemodynamic responses to orthostasis in normal subjects. *Aviat Space Environ Med* 2006;77(5):486-493.
18. **Cooke WH** et al. Human autonomic and cerebrovascular responses to inspiratory impedance. *J Trauma* 2006;60(6):1275-1283.
19. **Sigurdsson G** et al. Effects of an inspiratory ITD on blood pressure and short term survival in spontaneously breathing hypovolemic pigs. *Resuscitation* 2006 Mar;68(3):399-404.
20. **Smith SW** et al. Use of an ITD in hypotensive patients in the emergency department. *Circulation* 2006;114(Suppl II):18.
21. **Yannopoulos D** et al. Intrathoracic pressure regulation for intracranial pressure management in normovolemic and hypovolemic pigs. *Crit Care Med* 2006;34(12):S495-500.
22. **Jones AE** et al. Emergency department hypotension predicts sudden unexpected in-hospital mortality. *Chest* 2006;130:941-946.
23. **Convertino VA** et al. Inspiratory resistance maintains arterial pressure during central hypovolemia: implications for treatment of patients with severe hemorrhage. *Crit Care Med* 2007;35(4):1145-1152.
24. **Idris AH** et al. Imposed power of breathing associated with use of an ITD. *Respir Care* 2007;52(2):177-183.
25. **Melby DP** et al. Increased impedance to inspiration ameliorates hemodynamic changes associated with movement to upright posture in orthostatic hypotension: a randomized blinded pilot study. *Heart Rhythm* 2007;4(2):128-135.
26. **Metzger A** et al. An ITD improves 24-hour survival in a spontaneously breathing pediatric porcine model of hemorrhagic shock. *Circulation* 2007; Supplement II;116(16):II-632.
27. **Rickards CA** et al. Inspiratory resistance delays the reporting of symptoms with central hypovolemia: association with cerebral blood flow. *Am J Physiol Regul Integr Comp Physiol* 2007;293(1):R243-250.
28. **Eastridge BJ** et al. Hypotension begins at 110 mmHg: redefining "hypotension" with data. *J Trauma* 2007;63:291-299.
29. **Rickards CA** et al. Inspiratory resistance, cerebral blood flow velocity and symptoms of acute hypotension. *Av Sp Env Med* 2008;79(6):557-564.
30. **Ryan KL** et al. Breathing through an inspiratory threshold device improves stroke volume during central hypovolemia in humans. *J Appl Physiol* 2008;104:1402-1409.
31. **Voelckel WG** et al. Inspiratory ITD effects on hypotension in heat-stroked swine. *Av Sp Env Med* 2008;79(8):743-748.
32. **Parsons D** et al. The impedance threshold device (ITD-7): a new device for combat casualty care to augment circulation and blood pressure in hypotensive spontaneously breathing warfighters. *J Spec Ops Med* 2009;9(2):49-52.
33. **Metzger A** et al. Chapter 35: Harnessing cardiopulmonary interactions to improve circulation and outcomes after cardiac arrest and other states of low blood flow. *Handbook of Cardiac Anatomy, Physiology and Devices*. Humana Press 2009.
34. **Shih AC** et al. Cardiopulmonary effects of a new inspiratory impedance threshold device in anesthetized hypotensive dogs. *Vet Anaesth Analog* 2010;37:215-221.
35. **Vigani A** et al. Cardiopulmonary effects of a new inspiratory impedance threshold device in acute hemorrhagic shock in dogs. *J Vet Emerg Crit Care* 2011;21(6):618-624.
36. **Smith SW** et al. An ITD increases blood pressure in hypotensive patients. *J Emerg Med* 2011;41(5):549-558.
37. **Convertino VA** et al. Optimizing the respiratory pump: harnessing inspiratory resistance to treat systemic hypotension. *Respir Care* 2011;56(6):846-857.
38. **Suresh M** et al. Treatment of hypotensive patients by EMS personnel by modulating intrathoracic pressures. *Prehosp Emerg Care* 2012;16(1):173.
39. **Metzger A** et al. Augmentation of negative intrathoracic pressure improves hemodynamics without popping the clot: a randomized study comparing an impedance threshold device versus saline to treat severe hemorrhage through permissive hypotension in a spontaneously breathing porcine model. *Prehosp Emerg Care* 2012;16(1):174.
40. **Convertino VA** et al. Use of respiratory impedance in prehospital care of hypotensive patients associated with hemorrhage and trauma: a case series. *J Trauma Acute Care Surg* 2012; 73:S54-S59.
41. **Segal N** et al. Use of an impedance threshold device to treat severe hypotension in a pregnant woman: case report and review of the literature. *J Emerg Med* 2013; in press.

The generally cleared indication for the ResQGARD is for a temporary increase in blood circulation during emergency care, hospital, clinic and home use. Research is on-going in the United States to evaluate the long-term benefit of the ResQGARD for indications related to hypotension from specific etiologies. The studies listed here are not intended to imply specific outcome-based claims not yet cleared by the US Food and Drug Administration.



1905 County Road C West
Roseville, MN 55113 USA
Phone: 651-403-5600
Toll-free: 877-737-7763

info@advancedcirculatory.com
AdvancedCirculatory.com