

ALLEVAC BUCKLE RESTRAINT SYSTEM

\$54.60

SKU: 53-203

NSN: 6530-01-626-0965 |



PRODUCT DESCRIPTION

The Allevac Buckle Restraint System is a set of straps that provide chest, waist, thigh and ankle restraint. The unique Allevac Buckle is the superior restraint buckle. It allows for simple, effective operation and the Allevac Buckle cannot be broken or fouled by sand or debris like other buckles. The Allevac Buckle has been tested to a 2250 lb rating.

SPECIFICATIONS

- Tested to a 2250 lb rating
- Unit Dimensions 6' L x 3"W x .188"D
- Unit Cube .02
- Universal
- Color: Black