

**ALLEVAC® MANTIS SYSTEM****From: \$274.56**

Ultimate Evacuation System

SKU: 53-302**NSN:** 6530-01-620-2510 |**PRODUCT DESCRIPTION**

The AllEvac Mantis System is a compact, lightweight evacuation system containing the Tri-Fold Evacuation Backboard, Buckle Restraint Straps (5), Casualty Items Equipment Bag and AllEvac Mantis Carrier. It is the ultimate in simplicity and durability. This System is ideal for vehicle extraction due to its slimline design and lightweight configuration.

The Mantis Board is made of lightweight aluminum and has zero risk of fire hazards, cracking, warping and gouging that is common with plastic and composite boards in tactical and vehicle operations.

The AllEvac Mantis Carrier is made of 1000 denier Cordura nylon. It is Berry Compliant and available in four color options.

SPECIFICATIONS

- Specifications:
 - Dimensions: 26" x 17.5" x 5"
 - Unit Weight: 18 lbs
 - Weight Capacity: 500 lbs
- Carrier available in four color options
 - Coyote Brown, Black, Foliage Green, Multicam
- Berry Compliant
- Made in USA