

ALLEVAC® PATIENT RECOVERY BED**\$627.95**

Choose to be Organized

**SKU:** 57-000**NSN:** 6530-01-590-7741 |**PRODUCT DESCRIPTION****The Recovery Bed Solution for Mass Casualties**

The AllEvac® Patient Recovery Bed was designed to help provide a fast, effective, and economical response to mass casualty and hospital surge situations. The bed's lightweight, aluminum frame supports up to 450 lbs, enabling care to almost any size patient. In order to maximize patient comfort and treatment, this bed has a 5-position locking head elevator and 2-position foot elevator; this enables a modified Trendelenburg position. The AllEvac® Patient Recovery Bed Set comes complete with Safety Rails, IV Pole, Flat Bed Sheets (PK/2) and Safety & Instructional Labeling. It is intended for rapid deployment at fixed or portable locations with the ability to be stored securely and easily transported.

PRODUCT INFO

- 5-position locking head elevator enables modified Trendelenburg position
- Sets up in seconds
- Stores and transports easy with lightweight aluminum

SPECIFICATIONS

- Part Number: 57-000
- NSN: 6530-01-590-7741
- Specifications:
 - Working Capacity Weight: 450 lbs
 - Bed Weight: 30 lbs

- Deployed Dimensions: 80" x 32" x 18"
- Height adjustable at head and foot
- Wheelchair Transferrable Height (18")
- Antimicrobial Soil/Resistant Triple Layer Vinyl Mattress
- Spring and Deck Design (Reduces Bed Soars)
- Easy to Decontaminate
- ADA Compliant
- Made in USA