

**APLS LIFE GUARD****\$118.06**

Effective Portable Litter

SKU: 59-220**PRODUCT DESCRIPTION**

APLS Life Guard is an extremely effective portable litter that provides essential thermal and absorbent protection during emergency rescues and evacuations. APLS Life Guard provides the features of a more extensive portable litter system, but in a more efficient and streamlined design.

The compact portable litter features a Mylar lining that helps retain body heat while a treated nylon backing protects against wind and moisture. A highly absorbent cellulose core wicks fluids away from the patient, thus providing an additional level of protection. A robust Velcro® zipper provides easy access to the patient and eight padded handles enable first responders to quickly transport the patient.

PRODUCT INFO**Product Attributes**

- Built-in thermal properties for warmth during transport
- Treated nylon backing protects against wind and moisture
- Full-size litter is lightweight, rugged, strong
- Supports 320 - 350 lbs (145.15 - 158.76 Kg) wet or dry
- Detachable, Mylar-lined hood

Clinical Benefits

- Contains up to 1.25 liters of blood and bodily fluid
- Velcro zipper provides access to point of injury

- Eight reinforced handles enable easy carrying and transport
- Pocket for medical documents and personal effects

SPECIFICATIONS

- NSN: Pending
- PN: 59-220
- Not made with natural rubber latex / Single patient use
- Target Absorbency (H₂O): 1,000 - 1,450 grams, target 1,250 grams (1.25 liters)
- Pad Thickness: 0.15"(0.38 cm)
- Dimensions-
 - Deployed: 32" x 81" (81.28 cm x 205.74 cm)
 - Weight: 1.4 lbs. (.635 kg)
 - Weight Capacity: 320 - 350 lbs. (145.15 - 158.76 Kg)
- Made in USA