

**APLS RIGID THERMAL GUARD****\$632.37**

Rugged Litter with Thermal Protection

SKU: 59-225**NSN:** 6530-01-620-2667 |**PRODUCT DESCRIPTION**

APLS Rigid Thermal Guard is a lightweight portable litter that provides the structural benefits of a traditional military litter and improves the evacuation and transportation of injured personnel, particularly in remote, rugged locations. Weighing half of a traditional NATO litter, the expandable infrastructure and highly absorbent cellulose core form a truly unique, state-of-the-art litter system.

The absorbent core contains and controls blood and bodily fluids and a treated nylon backing protects against wind and moisture, thus increasing patient comfort while reducing the risk of cross contamination. Delivered as a compact 19" x 14" x 4" package, APLS Rigid Thermal Guard can be easily carried into the battlefield where it can expand quickly into a full-size rigid litter. For additional thermal protection, the APLS Rigid Thermal Guard features a Mylar lining that holds in heat and protects the patient against hypothermia during transport and care.

PRODUCT INFO**Capabilites**

- SUPPORT: 320 – 350 lbs. wet or dry
- THERMAL: Mylar lining for heat retention
- PORTABLE: litter includes rigid, expandable infrastructure system

Product Attributes

- Portable litter includes rigid, expandable infrastructure system
- Eight reinforced handles allow for easy carrying and transport

- Leak proof backing for additional strength while dragging
- Supports 320 - 350 lbs. (145.15 - 158.76 Kg) wet or dry

Clinical Benefits

- Contains up to 1.7 liters of blood and bodily fluids
- Mylar lining features thermal properties for heat retention
- Convenient close-able access points improve patient treatment
- Lightweight and rugged litter includes pocket for medical papers, personnel effects

SPECIFICATIONS

- NSN: 6530-01-620-2667
- PN: 59-225
- Not made with natural rubber latex / Single patient use
- Target Absorbency (H₂O): 1,500 - 1,850 grams (Target =1,700 grams, or 1.7 liters)
- Pad Thickness:0.7250" (1.84 cm)
- Dimensions-
 - Packaged: 19" x 14"x 4" (48.26 cm x 35.56 cm x 10.16 cm)
 - Weight: 10.0 lbs
 - Weight Capacity: 320 - 350 lbs. (145.15 - 158.76 Kg)
- Made in USA