

APLS SURGERY SOAKER/CATCH PAD**From: \$81.11**

Stay Dry and Clean



SKU: N/A**NSN:** 6530-01-561-7741 |**PRODUCT DESCRIPTION**

APLS Surgery Soaker is a disposable absorbent pad that contains blood, bodily fluids, and solutions under operating room tables in all military and civilian medical environments.

Constructed with a highly absorbent core and poly backing, APLS Surgery Soaker is the safe way to maintain dry and clean conditions in all medical environments.

APLS Surgery Soaker's fluid absorption properties enhance the safety, cleanliness and infection control related to patient surgery. The highly absorbent pad retains all solutions including saline, betadine, and sterilants.

APLS Surgery Soaker improves workplace safety while reducing the expense, time and risk of clean up.

PRODUCT INFO

- Contains up to 3.25 liters of blood and bodily fluid
- Leak proof poly backing
- Reinforced straps secure pad to operating room table
- Lightweight, strong and durable
- Supports 1.02 - 1.06 lbs. (0.46 - 0.48 kg)
- Compressed vacuum packed

SPECIFICATIONS

- Finished Product Size: 28" x 54" (71.12 cm x 137.16 cm)
- Target Absorbency (H₂O): 3,000 – 3,500 grams (Target 3,250 grams, or 3.25 liters)
- Pad Thickness: 0.0725" (0.18 cm)
- Pad Weight: 2.74 lbs. (1.24 kg)