

APLS THERMAL GUARD SHEATH**From: \$328.90**

SKU: N/A**PRODUCT DESCRIPTION**

APLS Thermal Guard Sheath is designed to work in concert with standard military litters and improves patient warmth during transient transport. A fleece-lined panel protects the patient against hypothermia while a treated nylon backing protects against wind and moisture.

Symmetric, fleece-lined hood openings at each end of the litter allow the patient to be placed in the APLS Thermal Guard Sheath from either side.

The highly absorbent core contains blood and bodily fluids, thus reducing the risk of cross-contamination.

APLS Thermal Guard Sheath can be used independently, on top of, or as an envelope around a structured litter.

Each APLS pad includes eight reinforced handles that allow for easy carrying and transport.

PRODUCT INFO

- Access points on both sides improve patient treatment
- Symmetric design enables access from both ends
- Built-in thermal properties for warmth during transport
- Contains up to 4.5 liters of blood and bodily fluids
- 8 padded handles allows for easy carrying and transport
- Supports 420 - 450 lbs. (190.51 - 204.12 kg) wet or dry
- 4 clear pockets on top for storage of treatment products

- Leak proof nylon backing
- Lightweight but strong
- Compressed vacuum packed

SPECIFICATIONS

- **Finished Product Size:** 33" x 81" (83.82 cm x 205.74 cm)
- **Target Absorbency (H₂O):** 4,250 - 4,750 grams (target = 4,500 grams, or 4.5 liters)
- **Pad Thickness:** 0.3150" (0.80 cm)
- **Pad Weight:** 8.20 - 8.80 lbs. (3.72 - 3.99 kg)