

**APLS THERMAL GUARD****\$406.51**

Highly Absorbent

SKU: 59-210**NSN:** 6530-01-562-6008 |**PRODUCT DESCRIPTION**

APLS Thermal Guard is a portable and disposable absorbent litter that provides maximum thermal protection during emergency transport. APLS Thermal Guard provides an integrated system for improved patient care: the highly absorbent cellulose core wicks fluids away from the patient, the fleece-lined bag and detachable hood retain body heat, and the treated nylon backing protects against wind and moisture.

Two close-able “L” shaped access points enable the patient to retain body heat while being treated by medical personnel. For added thermal protection, optional internal pockets provide a “heat zone” in four key areas when used with a mobile heat source.

APLS Thermal Guard can retain 4.5 liters of bodily fluid while safely supporting up to 450 pounds. Two pockets inside the bag hold patient documents, x-ray disks and personal effects.

PRODUCT INFO**Capabilities**

- **SUPPORT:** 420 – 450 lbs. (190.51 -204.12 kg) wet or dry
- **FLEXIBLE:** Detachable hood components
- **DURABLE:** Lightweight and rugged design

Product Attributes

- Built-in thermal properties for warmth during transport
- Contains up to 4.5 liters of blood and bodily fluid

- Leak proof nylon backing
- Closable access points aide with patient treatment and warmth

Clinical Benefits

- Lightweight, rugged, strong
- Eight handles for easy transport
- Bottom straps secure to litters
- Locker for patient documents
- Detachable hood component

SPECIFICATIONS

- NSN: 6530-01-562-6008 (case of 5)
- PN: 59-210
- Not made with natural rubber latex / Single patient use
- Target Absorbency (H₂O):4,500 to 5,000 grams (Target= 4,750 grams, or 4.5 liters)
- Pad Thickness: 0.2250" (0.57 cm)
- Dimensions-
 - Packaged: 33" x 91"(83.82 cm x 231.14 cm)
 - Weight: 6.60 - 6.65 lbs. dryweight (2.99 - 3.2 kg)
 - Weight Capacity: 420 - 450 lbs. (190.51 -204.12 kg) wet or dry