



## **APLS TRANSPORT PAD WITH COMFORT PLUS**

**From: \$102.72**

Comfort is a plus

---

**SKU:** N/A

**NSN:** 6530-01-598-4810 |

### **PRODUCT DESCRIPTION**

APLS Transport Pad with Comfort Plus is a portable and disposable litter that includes all of the features of our signature APLS Transport stretcher plus an additional layer of padding for incremental patient comfort. APLS Transport with Comfort Plus is constructed with a durable canvas and highly absorbent core that absorbs up to 2.5 liters of fluid and supports up to 325 pounds, wet or dry.

The additional padding provides comfort to pressure points while the treated nylon backing protects against wind and moisture.

APLS Transport with Comfort Plus can be used for emergency patient transport, on top of hard litters and stretchers, in operating rooms and as a transfer sheet.

APLS Transport with Comfort Plus is available with an anti-microbial treatment that helps control odors caused by body fluids.

### **PRODUCT INFO**

- Foam-like padding contained inside litter
- Available with anti-microbial treatment for odor control
- Lightweight, rugged and portable
- Unique absorbent core contains up to 2.5 liters of blood and bodily fluids
- Built-in thermal properties for warmth during transport
- Eight (8) handles allow for carrying and easy transport
- Leak proof nylon backing suitable for dragging
- Straps to ensure patient stability

- Supports 275 - 325 lbs. (124.74 - 147.42 kg) wet or dry
- Compressed vacuum packed or bulk packed

## **SPECIFICATIONS**

- **Finished Size:** 20.0" x 72.0" (50.80 cm x 182.88 cm)
- **Target Absorbency (H<sub>2</sub>O):** 2,350 - 2,650 grams (Target = 2,500 grams, or 2.5 liters)
- **Pad Thickness:** 0.75" (1.91 cm) - Uncompressed
- **Pad Weight:** 2.81 lbs. (1.27 Kg)