

**BLOODBOXX™ BASIC NESTING SYSTEM****\$988.90**

FWB Transport System

SKU: 39-450**PRODUCT DESCRIPTION**

The BloodBoxx Basic Nesting System is a modular, fresh whole blood (FWB) transport system that provides thermal protection for up to 6 bulk FWB units or 2 nested BloodBoxx Tactical units with no power requirements. The system is designed to enable safe and cost-effective transport of FWB for 72-96 hours* while maintaining the required temperature range of 1°C to 10°C .

By combining state-of-the art, ultra efficient insulated chambers with optimized thermal management materials, the BloodBoxx Basic Nesting System delivers the controlled temperature transport solution to better support unique missions in austere environments.

*Temperature maintenance time varies based on payload configuration (Nested units vs. Bulk FWB units)

PRODUCT INFO

- 72-96 hours of transport time
- The modular nesting system allows you to carry FWB to an austere, forward staging area and conduct follow-on disbursed operations.
- Maintains FWB at proper temperature if there is an unexpected delay in mission execution.
- Modular payload configurations
- No power source required
- Carrying FWB means no immediate, if any, need to initiate a walking blood bank.
- Minimizes individual fighting load while increasing distributed capability on the target.

- Designed for 500mL CDPA-1 Whole Blood Bags

SPECIFICATIONS

- PN: 39-450 (System) – Modular Payloads at 1-10C
 - 4-6 x 500ml FWB for 72-96 hours
 - 2 x Single Unit Boxes (Tacticals) for 72-96 hours
 - Dimensions: 12.5”L x 13.5”W x 11”H
 - System Weight: 14.8 lbs **
 - Weight with 2 x Tacticals: 17.2 lbs **
- PN: 39-101 (Tactical)
 - 1 x 500ml FWB
 - Maintains 1-10C for 12-24 hours of mission time
 - Dimensions: 5”L x 3 3/4”W x 7 3/8”H
 - Weight: 1.7 lbs **
- 5 Year Shelf Life
- Berry Compliant
- Reusable
- Made in USA

** Listed weights do not include FWB payload