

**BLOODBOXX™ TACTICAL****\$324.50**

Blood Carried Far Forward

SKU: 39-101**PRODUCT DESCRIPTION**

Our enhanced temperature control technology delivers the next evolution for thermo-regulated payload requirements. By combining state of the art, ultra-efficient, insulated chambers with optimized thermal management materials, we offer the controlled temperature transport solution to better support unique missions in austere environments. The BloodBoxx Tactical delivers a single unit of blood to the battlefield in a durable, tested system. Using no additional power, this revolutionary design maintains the required temperature range for up to 24 hours, depending on ambient temperature.

PRODUCT INFO

- WEIGHT: Less than 2lbs
- DURABLE: 1000 denier shell
- TIME: 12-24 hours of transport

Product Attributes

- Light weight
- Low Cube
- No power requirements
- Simple reusable system
- Optimization of thermal management materials
- Maintains between 1°C and 10°C
- Designed for 500ml CDPA-1 Whole Blood Bags

SPECIFICATIONS

- PN: 39-101
- ISTA Standard Tested
- 5-year Shelf Life
- Specification:
 - Unit Dimensions: 5"L x 3 3/4"W x 7 3/8"H
 - Unit Weight: 1.7 lbs
- Available Direct
- Berry Compliant
- Made in USA