

CELOX™ RAPID RIBBON

From: \$23.05

Narrow track wound treatment

SKU: 30-223

NSN: 6510-01-623-9442 |



PRODUCT DESCRIPTION

Celox™ Rapid Ribbon brings all the benefits of Celox Rapid Gauze in a 1-inch strip. This Chitosan impregnated dressing is designed to stop moderate-to-severe bleeding with one minute of compression time. Its low-profile packaging enables Celox Rapid Ribbon to be easily integrated into existing First Aid Kits. Its increased capability is ideal for the treatment of narrow track wounds such as bullet and stab injuries. Celox Rapid Ribbon follows the MARCH mnemonic in a color-coded, trauma treatment system distinctly packaged with icon identifiers to treat the five leading causes of preventable death following traumatic injury.

PRODUCT INFO

- 2 Minutes Faster than Closest Competitor
- Cleared by FDA for 1 Minute of Compression
- Thin strip configuration
- Reduces Blood Loss
- Works in hypothermic conditions
- Treatment of narrow track wounds
- Lightweight
- Low Cube

SPECIFICATIONS

- PN: 30-223 (EA)
 - NSN: 6510-01-623-9442
- PN: 30-223-150 (CS/150)

- NSN: 6510-01-623-9730
- Sterile, Single-Patient Use
- 5-year shelf life
- Dimensions (EA):
 - Packaged: 5.25"L x 3.5"H x 0.5" W
 - Useable: 1"W x 5'L
- Weight (EA): 0.5 oz
- TAA Compliant