

COMPRESSED ROLLED GAUZE

From: \$3.86

SKU: N/A



PRODUCT DESCRIPTION

Our Compressed Rolled Gauze takes care of moderate to severe wounds involving heavy blood loss. Made with high-quality cotton for maximum absorption and stability, Compressed Rolled Gauze is vacuum packed small for easy storage yet expands to cover large wound areas. The product comes with a 4.5" x 4.1 yards (stretched) roll of crinkle fluff USP Type VII compliant gauze, sterile and ready to use.

Compressed Rolled Gauze is the perfect companion to our compression bandages to stop severe blood loss and save lives.

PRODUCT INFO

- Item Number: HHPG or HHPG-CA
- Quantity: Individual or CS/175

- **HHPG NSN: 6510-01-503-2117**

SPECIFICATIONS

- 6-Ply Cotton
- Sterile
- Primed Gauze
- TAA Compliant
- Packaged Dimensions: 3"x2"x1"
- Deployed Dimensions: 4.5"x4.1 yards stretched
- Weight: 2 oz
- Shelf Life: 7 Years