

**CRAVAT****\$2.65**

Versatile Textile - Dependable Treatment

SKU: 32-600**NSN:** 6510-01-669-7885 |**PRODUCT DESCRIPTION**

The Cravat is a multipurpose triangular shaped textile that can be used to sling and support injuries, immobilize fractures and dislocations, or even bandage wounds. The bandage is packaged in a low cube configuration and is made of Berry Compliant muslin cloth. The Cravat easily fits into any first aid kit and is a staple for basic bandaging or splinting.

PRODUCT INFO**Capabilities**

- **Versatile** Treatment Capability
- **Secures** and Stabilizes
- **Simple** Instructions for Use

Product Attributes

- 37" x 37" x 52" Triangular Bandage
- 2 Safety Pins to Secure Placement
- MARCH™ product labeling and packaging

Clinical Benefits

- Product Uses:
 - Sling
 - Immobilization / Splint Support
 - Medical Wrap / Bandage

SPECIFICATIONS

- NSN: 6510-01-669-7885
- PN: 32-600
- Specifications:
 - Unit Dimensions:
 - Packaged: 3.5" x 2" x .625"
 - Deployed: 37" x 37" x 52"
 - Unit Weight: 1.6 oz
- Available Direct, Prime Vendor, ECAT and GSA
- Berry Compliant
- Made in USA