

**DART™****\$10.11**

The Pre-Hospital Choice

**SKU:** 21-001**NSN:** 6515-01-601-1292 |**PRODUCT DESCRIPTION**

The MARCH™ Dart™ gives providers at any skill level the confidence and capability to perform reliable chest decompression and improve outcomes for casualties with a tension pneumothorax, the third leading cause of preventable traumatic death.<sup>1</sup>

The Dart™ is a sterile, 14 gauge x 3.25" needle with a uniformed, thick walled catheter. The design minimizes risk of catheter collapse in the intercostal space. The sure-grip hub ensures stability in difficult environments. The rugged case was designed for easy, one-handed opening and doubles as a sharps container for patient and practitioner safety.

The Dart's™ six-year shelf life lowers lifecycle costs when compared to most chest decompression needles.

**PRODUCT INFO****Capabilities**

- Designed with rugged case, sure-grip hub, and thick-walled catheter
- Minimal Risk for additional decompression attempts
- Improve Outcomes when treating tension pneumothorax

**Product Attributes**

- Ergonomic sure grip hub:
  - Easy to use in difficult environments

- Easy grip with gloves
- Snug fit in case – no falling out
- Thick-walled, uniform catheter
- Patented sterilization plug
- Case doubles as a sharps container
- Fits in a military IFAK
- Low volume, low cube

## Clinical Benefits

- Improves outcomes when treating tension pneumothorax, the third leading cause of preventable battlefield death<sup>1</sup>
- Catheter design minimizes risk of clogging, kinking and need for additional decompression attempts
- Easy access case for quick use and safe needle disposal
- Purpose built decompression needle, not a re-purposed IV Catheter

## SPECIFICATIONS

- NSN: 6515-01-601-1292
- PN: 21-001
- FDA Listed, ISO 13485 and cGMP Compliant
- Sterile / Single-use
- 6-year shelf life
- Specifications:
  - Unit Dimensions
    - 6.25”H x 0.88”W x 0.5”D
    - Needle Dimensions: 14 gauge x 3.25”
  - Unit Weight: 0.5 oz
- Not made with natural rubber latex
- Available Direct, Prime Vendor, ECAT, CEC and GSA
- 100% TAA Compliant
- Made in USA

## TRAINING

[Dart IFU](#)

## REFERENCES

**1** Eastridge BJ, et al., Death on the battlefield (2001-2011): implications for the future of combat casualty care., J Trauma Acute Care Surg. 2012 Dec;73(6 Suppl 5):S431-7.