



ELDONCARD® BLOOD TYPING KIT

\$8.25

Know Your Blood Type

SKU: 80-810

NSN: 6550-01-587-1889 |

PRODUCT DESCRIPTION

The EldonCard® Blood Typing Kit contains everything required to determine blood group and Rh factor using a whole blood sample obtained from a finger prick without the need for separate reagents or laboratory training. The EldonCard is a patented test card device that comes specially treated and impregnated with dried antibody sera Anti-A, Anti-B and Anti-Rh used universally for determining blood type by the forward typing method. This kit creates a permanent Blood Type Record Card for medical histories and emergencies.

The kit contains: ELDONCARD 2511, Alcohol preparation swab, EldonSticks (4), Auto-Lancet (2), ELDON FOIL, Plastic dropper and Instructions For Use.

SPECIFICATIONS

- NSN: 6550-01-587-1889
- Part Number: 80-810
- Unit of Issue: Each
- CAGE Code: 515E6
- 24 month shelf life
- Single Use
- Results after just 2 minutes
- Dimensions: 6" x 4.25" x 0.5"
- 100% TAA Compliant