



## **EMMA® CAPNOGRAPH**

**From: \$250.00**

---

**SKU:** 25-600

**Info::** Airway Adapter sold separately |

## **PRODUCT DESCRIPTION**

EMMA® is a portable, real-time capnography device that measures end-tidal carbon dioxide and respiration rate. Designed to fit easily onto a breathing circuit, EMMA can be utilized at multiple points of care including pre-hospital, emergency medicine, operating room, intensive unit, and long-term acute care. Its components are easy to maintain with no calibration required and it has a rugged design that makes it reliable for operation in challenging environments.

## **PRODUCT INFO**

- Part Number: 25-600
- Dimensions: 2.1 x 1.5 x 1.5"
- Weight: 2.1 oz.
- Includes EMMA, pouch and lanyard
- Battery Life Up to 10 Hours
- Audible and Visual Alarm System
- Minimal warm-up time
- No Calibration Required
- 14.4 Second Sweep of CO2 Values
- Not made with natural rubber latex
- Made in the USA