

EZ-IO® MANUAL NEEDLE SET**\$1,500.00**

SKU: 35-094**Info::** Quantity: Box of 6 |**PRODUCT DESCRIPTION**

The Arrow® EZ-IO® Manual Needle Set is designed to provide controlled vascular access via the intraosseous (IO) route to patients in emergent situations when standard venous access is not possible. Using IO access, providers can deliver lifesaving fluids and medications rapidly and effectively.^{1 2}

PRODUCT INFO

- New longer, 38.5 mm needle (previously 25 mm)
- Indicated for six insertion sites (proximal humerus and proximal and distal tibia, bilaterally)
- Provides controlled IO access in emergent situations when standard venous access is not possible
- Intuitively designed for pre-hospital providers
- Manual insertion - no driver (or powered driver) required
- Easily stored in a medic bag
- Set includes: EZ-IO® Manual Needle, EZ-Stabilizer® Dressing, EZ-Connect® Extension Set, NeedleWISE® Sharps Block & Wristband, and Quick Reference Guide

SPECIFICATIONS

- PN: 35-094
- FDA Cleared
- Specifications:
 - Dimensions: 8.3" x 5.1" x 2.0"
 - Weight: 60.7g

- 4 year shelf life
- Sterile
- Latex-Free
- 304 Stainless Steel Needle
- 100% Made in USA

REFERENCES

1. Lewis P, Wright C. Saving the critically injured trauma patient: a retrospective analysis of 1000 uses of intraosseous access. *Emerg Med J* 2014;31(9):784. doi:10.1136/emered-2014-203588
2. Paxton JH, Knuth TE, Klausner HA. Proximal humerus intraosseous infusion: a preferred emergency venous access. *J Trauma* 2009; 67: 606-11. Research sponsored by Teleflex.