

EZ-IO® MANUAL TRAINING NEEDLE SET

\$1,155.00

*Image depicts the non-training EZ-IO Manual Needle. EZ-IO Manual Training Needle comes with **red colored** hub.



SKU: 35-094T

Info:: Quantity: Box of 6 |

PRODUCT DESCRIPTION

The EZ-IO® Manual Needle Training Set allows medical personnel to train in six (6) FDA-cleared peripheral sites for the immediate manual securing of vascular access. Its unique, red colored hub provides an immediate visual cue that it is only intended for training purposes. The manual trainer needle version features the same weight and dimensions as the tactical use version so trainees and medical personnel can familiarize themselves with the proper use and techniques required for in-field applications.

PRODUCT INFO

- New longer, 38.5 mm needle (previously 25 mm)
- Enables training in six insertion sites (proximal humerus and proximal and distal tibia, bilaterally)
- Unique, red colored hub indicates for training use only
- Manual insertion-no driver (or powered driver) required
- Set includes: EZ-IO® Manual Training Needle, EZ-Stabilizer® Dressing, EZ-Connect® Extension Set, NeedleVISE® Sharps Block & Wristband, and Quick Reference Guide

SPECIFICATIONS

- PN: 35-094T
- FDA Cleared
- Specifications:
 - Dimensions: 8.3" x 5.1" x 2.0"
 - Weight: 60.7g

- Non-Sterile
- Not For Human Use
- Latex-Free
- 304 Stainless Steel Needle
- 100% Made in USA