

EZ-IO® TRAINING NEEDLE SETS

From: \$51.25

EZ-IO Training Needles

SKU: N/A



PRODUCT DESCRIPTION

The Arrow® EZ-IO® Intraosseous Vascular Access System creates a stable port into the intraosseous space to access the central vascular system within seconds, causing minimal trauma to the patient. These EZ-IO Training Needle Sets are intended for use with the EZ-IO Power Driver.

SPECIFICATIONS

These Training Needle Sets are intended for use with the Arrow® EZ-IO® Intraosseous Vascular Access System Power Driver. The 15-gauge IO needles are available in 15mm (pediatric), 25mm (adult) and 45 mm (large adult) lengths. The 25mm needles are specially designed with a beveled drill tip to quickly and easily penetrate through the hard exterior of adult bones, while the 15mm needle permits easy access through the softer bones of pediatric patients. The new 45mm needle is designed for patients whom have excess tissue over the insertion site.

EZ-IO 25mm (Adult) Training Needle Set:

NSN: 6515-01-540-9808

- (1) 25mm Training IO Needle
- (1) Directions
- (1) EZ-Connect
- (1) Wristband

EZ-IO 15mm (Pediatric) Training Needle Set

NSN: 6515-01-541-9382

- (1) 15mm Training IO Needle
- (1) Directions

- (1) EZ-Connect
- (1) Wristband

EZ-IO 45mm (Large Adult) Training Needle Set:

NSN: 6510-01-597-0531

- (1) 45mm Training IO Needle
- (1) Directions
- (1) EZ-Connect
- (1) Wristband