



FRESH WHOLE BLOOD TRANSFUSION SET

\$106.70

Fresh Whole Blood Anywhere

SKU: 80-801

NSN: 6515-01-657-4750 |

PRODUCT DESCRIPTION

The Fresh Whole Blood Transfusion Set contains everything needed to test for ABO/Rh compatibility and perform a fresh whole blood transfusion in a compact sealed package and was developed with extensive research and feedback from the Armed Services Blood Program Office.

The transfusion of fresh whole blood is the best option for management of casualties with hemorrhagic shock. Blood transfusions have been used successfully in military conflicts and other austere conditions to resuscitate trauma patients suffering from shock due to hemorrhage.

This self contained set is an affordable and sensible solution for training and equipping organizations required to collect and transfuse fresh whole blood.

NOTE: Only personnel trained in blood collection, transfusion, and management of transfusion reactions should use this set.

PRODUCT INFO

Product Attributes

- 100% TAA Compliant Components
- Single Lot Tracking / UDI Compliant
- Self Contained Donor and Recipient specific modules in tamper evident re-sealable bags
- 12 month self life

Clinical Benefits

- Adhesive blood bag label for recording donor information
- EldonCard™: ABO/Rh screening card for simple, accurate, and fast results
- Re-sealable bag used to send all supplies back with casualty for retrospective testing
- Assembled in FDA Registered Facility

SPECIFICATIONS

- Part Number 80-801
- NSN 6515-01-657-4750
- CAGE Code: 515E6
- FDA Listed / UDI Compliant (Class II Device)
- 12 month shelf life
- Exceeds ASBP CPGs for Walking Blood Banks
- Specifications:
 - Package: 10.5"W x 14.0"H x 2.0"D
 - Weight: 18.5 oz
- 100% TAA Compliant
- Available for purchase: Direct, PV and ECAT

REFERENCES

- Part Number 80-801
- NSN 6515-01-657-4750
- CAGE Code: 515E6
- FDA Listed / UDI Compliant (Class II Device)
- 12 month shelf life
- Exceeds ASBP CPGs for Walking Blood Banks
- Specifications:
 - Package: 10.5"W x 14.0"H x 2.0"D
 - Weight: 18.5 oz
- 100% TAA Compliant
- Available for purchase: Direct, PV and ECAT