

**G3 POWER TRAINING DRIVER****\$242.00**

SKU: 35-059**NSN:** 6515-01-596-8614 |**PRODUCT DESCRIPTION**

The G3 Training Power Driver is in a fully sealed external case that houses the motor and Lithium Batteries. Provide 10 years of maintenance free shelf life and a minimum of 500 insertions. Trigger guard to protect against inadvertent activation of the driver when stored and carried in tightly packed carts, bags or combat medical gear. The trigger guard is attached to the Power Driver with a lanyard eliminating the need for a specialized case. G3 Power Driver enhancements include a battery indicator that provides a visual warning for the operator. Designed to facilitate rapid readiness testing, the LED is non-flashing green when the batteries are at optimal performance and flashing red when 10% of the life battery life remains. Unit Dimensions: 6.5"L x 4.5"W x 2.5"D Unit Cube: .04