



HAWK™ ADVANCED HYPOTHERMIA MANAGEMENT SET

\$159.02

Hypothermia has met its match

SKU: 59-320

NSN: 6530-01-673-4063 |

PRODUCT DESCRIPTION

The *HAWK™* Advanced Hypothermia Management Set is designed to help prevent and manage hypothermia by combining the advanced technology of adjustable, air-activated warming grids and the Thermal Patient Litter System.

The APLS Life Guard Litter System is a portable, disposable litter that provides essential thermal protection during emergency transport. The compact litter features a Mylar lining that helps retain body heat while a highly absorbent pad wicks fluids away from the patient, providing an additional level of protection.

The *HAWK™* Warming Grid safely activates in just minutes, using advanced technology to maintain ideal temperature. Its simple shut-off capability also extends the life of the warming grid while enabling temperature control.

PRODUCT INFO

- Produces constant heat 4-32 hours, depending on therapeutic needs
- Advanced technology maintains ideal temperature
- Safe, lightweight, low profile
- Quick, air-activated, adjustable heating process
- Easy on-off capability
- Thermal Patient Litter System supports 320 lbs (wet) – 350 lbs (dry)
- Built-in thermal properties for warmth during transport
- Vacuum sealed, ruggedized package
- Treated nylon backing protects against wind and moisture

- Detachable Mylar-lined hood
- Retains up to 1.5 liters of blood and bodily fluid

SPECIFICATIONS

- PN: 59-320
- NSN: 6530-01-673-4063
- Single Patient Use
- Not made with natural rubber latex
- Target Absorbency (H₂O): 1,000 – 1,450 grams, target 1,250 grams (1.25 liters)
- Pad Thickness: 0.15" (0.38 cm)
- Specifications:
 - Package Dimensions: 9.5"L x 7"H x 3.5"W
 - Weight: 2.10 lbs
 - Weight Capacity: 320 lbs (wet) – 350 lbs (dry)
- Made in USA