



IV INFUSION STARTER SET

\$24.99

Very Fast Venous Access

SKU: 80-100

NSN: 6515-01-633-0975 |

PRODUCT DESCRIPTION

The IV Infusion Starter Set allows for quick and effective extremity venous cannulation under the most demanding conditions in one complete set. It contains the necessary items for starting and securing a saline lock plus a macro-drip IV administration set. The rugged, low-cube, vacuum-sealed pouch keeps the contents of the kit in prime condition until needed.

PRODUCT INFO

IV Infusion Starter Set

Product Number:	80-100
Components	QTY
18 GA Catheter	1
20 GA Catheter	1
Saline Flush (5cc pre-filled syringe)	1
Povidone Pad	1
Alcohol Prep Pad	2
IV Administration Set	1
Needle-Free Valve w/ Luer Lock	1
Constricting Band	1
IV Catheter Dressing	1
Transparent Film Dressing	1
Adhesive Bandage	1
2"x2" Gauze Pad (2 pack)	1

SPECIFICATIONS

- NSN: 6515-01-633-0975
- PN: 80-100
- Shelf Life: 6 - 18 months
- Single Patient Use
- Not made with natural rubber latex
- Specifications:
 - Dimensions: 3.75"W x 7"H x 1"D

- Weight: 2.8 oz
- Available: ECAT, PV, VIPA, GSA, CEC
- Made in USA