



## MIGHTYSAT™ RX FINGERTIP PULSE OXIMETER

**\$355.00**

---

**SKU: 42-700**

### PRODUCT DESCRIPTION

The MightySat™ Rx Fingertip Pulse Oximeter measures SpO2, PR, Pi, and PVi with its at-a-glance high resolution display. Its flexible finger pad enables secure application and minimizes external interference. The touch pad has a high resolution display with a rotating and adjustable color screen for clear display of all measurements. The MightySat™ Rx offers a unique solution to point of care treatment and initial patient diagnostics.

### PRODUCT INFO

- High Resolution Display
- Signal I.Q.® provides measurement confidence of SpO2
- Real time graphical display of patient's Pleth Waveform
- Touchpad Allows for Easy Access of Settings and Features
- Optional Bluetooth LE Capability
- Flexible Finger Pad
- Durable and Lightweight Design

### SPECIFICATIONS

- PN: 42-700
- Specifications:
  - Weight: .16 lbs.
  - 2.9" x 1.6" x 1.2"
- Operates with 2 (1.5) volt AAA batteries
- Battery life of approx. 1800 spot checks

- Operating Temperature: 41-104F
- Storage Temperature: -40 to 158F
- Operating Humidity 10 to 95%, non-condensing