

MOJO® B-CON BELT KIT**From: \$82.00**

SKU: N/A**PRODUCT DESCRIPTION**

The B-Con Belt Kit is a tactically stocked bleeding control kit. This compact, lightweight design attaches to a standard 2" belt and is easily accessible during stressful situations. The B-Con Belt Kit is intended for initial responders as a means of immediate treatment following traumatic injury. It uses MARCH™ medical gear in a color-coded trauma treatment system that treats the five leading causes of preventable death. This kit is offered in two different configurations to accommodate the needs of any department or situation.

PRODUCT INFO**SET COMPONENTS:****Mojo B-Con Belt Kit - Basic (PN 87-531-RC)**

- Celox Rapid B-Con
- TMT Tourniquet
- Pr. Exam Gloves
- Mojo Nylon Belt Pouch (Red, check pattern)
- MARCH Quick Reference Cards

Mojo B-Con Belt Kit - Advanced (PN 87-535-RC)

- Celox Rapid
- TMT Tourniquet
- Pr. Exam Gloves

- Mojo Nylon Belt Pouch (Red, check pattern)
- MARCH Quick Reference Cards

SPECIFICATIONS

- Weight: 0.35 lbs.
- Dimensions: 6" L X 3.5" W X 2.5" H
- Red check pattern
- Up to 5 year shelf life
- Made in the USA