



MOJO® BELT TRAUMA KIT

From: \$126.45

Small Design, Big Capabilities

SKU: 88-104

PRODUCT DESCRIPTION

Mojo® Belt Trauma Kit is a comprehensive, compact, trauma treatment kit configured in a rugged, nylon pouch. It's low-profile design and convenient zipper pouch make this kit ideal for casualty care following traumatic injury. It attaches to standard 2" utility and gun belts for easy access. The Belt Trauma Kit contains MARCH™ medical gear in a color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five leading causes of traumatic death. This kit is offered in three different content configurations to accommodate the needs of any department.

PRODUCT INFO








- Treats the Leading Causes of Traumatic Death
- Ergonomically Designed for Ease of Use
- 5 Year Shelf Life

SPECIFICATIONS

Specifications

- **PN: 88-104-01** (Basic)
- **PN: 88-104-02** (Intermediate)
- **PN: 88-104-03** (Advanced)
- Specifications:
 - Weight: 11.8 oz
 - Dimensions: 6.5" x 4.5" x 2"

- Up to 5 Year Shelf Life
- Single-Use
- Berry Compliant Nylon Pouch
- 100% TAA Compliant

MARCH™		Belt Trauma Kit	Belt Trauma Kit	Belt Trauma Kit
Color-coded trauma treatment system		Basic	Intermediate	Advanced
	Product Number:	88-104-01	88-104-02	88-104-03
	Components	QTY	QTY	QTY
	Celox Rapid	1	1	1
	TMT Tourniquet		1	1
	SWAT-T Tourniquet	1		
	Positional Airway Card	1	1	1
	FoxSeal	1	1	1
	Battle Wrap			1
	Emergency Blanket	1	1	1
	Belt Trauma Pouch	1	1	1
	Pair Exam Gloves	1	1	1
	Permanent Marker	1	1	1
	MARCH QuickStart Guide	1	1	1