

MOJO® COMBAT LIFESAVER BAG**From: \$188.67**

Bring Your Best as CLS

SKU: 83-320**PRODUCT DESCRIPTION**

The Mojo® Combat Lifesaver Bag uses MARCH™ medical gear in a user friendly, color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five preventable causes of death on the battlefield. The Mojo® CLS bag is mission designed and highly versatile. It contains components for a warfighter with advanced individual or CLS training to treat multiple casualties at point of injury. The Mojo CLS has three methods of attachment: adjustable shoulder strap, attach directly to MOLLE, or worn with the Mojo® MOLLE belt. The bag has a top zipper panel for full view of and quick access to its contents. It is made of 1000D Cordura and is available in four different color options (coyote brown, black, foliage green, multicam).

- Basic- PN 83-324-01
- Intermediate- PN 83-323
- Advanced- PN 83-325

Image shown is PN 83-325, actual product may vary

The Mojo® Combat Lifesaver Bag (CLS) Mirrors the USAMMA Approved 2009 CLS Packing List**PRODUCT INFO**

- Quick Access
- Intuitive Layout

- Versatile Portability

Product Attributes








- Ergonomic design evenly distributes weight and is worn comfortably on your body during long missions
- Wearable in three ways (belt, shoulder strap, MOLLE system)

Clinical Benefits

- Systematic access and full visualization of your life saving equipment ensures the smoothest casualty intervention
- 180-degree top-opening zipper panel for quick grab access
- Fully mission-configurable kit

SPECIFICATIONS

- PN: 83-325
- coyote brown, black, foliage green, multicam
- Up to 5-Year Shelf Life
- Single Use
- Dimensions: 16" x 9.5" x 6"
- Weight: 8.9 lbs
- Berry Amendment / Trade Compliant
- Made of 1000D Cordura
- Made In USA

MARCH™ Color-coded trauma treatment system		Mojo CLS Advanced	Mojo CLS Intermediate	Mojo CLS Basic
	Product Number:	83-325	83-323	83-324-01
	Components	QTY	QTY	QTY
	Celox Rapid	2	2	2
	TMT Tourniquet	2	2	2
	ITClamp	2		
	Prelubricated Nasal Airway	2	1	1
	Sentinel Chest Seal	2	2	2
	Dart			2
	Chest Decompression Set	2	2	
	Compression Dressing			2
	Elastic Bandage 6" x 4.5 yd			1
	Abdominal Dressing	1	1	1
	Battle Bandage	2	2	
	Battle Wrap	2	1	
	Rolled Gauze	2	2	2
	Polycarbonate Eye Shield w/ Garter	2		
	Aluminum Eye Shield w/ Garter		1	1
	Blizzard Light Blanket		1	1
	HAWK Warming Grid	1	1	
	APLS LifeGuard Military	1		
	Mojo CLS Bag	1	1	1
	Mojo Cravat	3	3	3
	Hook Knife	1	1	
	Pair Exam Gloves	4	4	4
	Universal Splint	1	1	1
	Cloth Medical Tape 3" x 10 yd	1	1	1
	Trauma Shears	1	1	1
	Shears Retractor	2	2	
	TCCC Card DD 1380	2	2	2
	Permanent Marker	2	2	2
	MARCH QuickStart Guide	1	1	1

NSNS

NSNs for Bag only:

Black: NSN 8465016657999

Coyote Brown: NSN 8465016658011