








MOJO® CRISIS RESPONDER**From: \$100.10**

SKU: N/A**PRODUCT DESCRIPTION**

The Crisis Responder is a tactically stocked trauma kit for First Responders. Its small-scale and compact design make it easy to access products during stressful situations. The Crisis Responder is intended for first responders and rescue task force members in emergency situations. It uses MARCH™ medical gear with icon-identifiers. This color-coded trauma treatment system addresses the five leading causes of preventable death following trauma. The Crisis Responder is offered in two different configurations to accommodate the needs of any department or situation.

SPECIFICATIONS

- Weight: 1.95 - 2.35 lbs
- Dimensions: 10" L X 7.25" W X 5"H
- Available in red or black diamond, check pattern
- Up to 5 year shelf life
- Made in the USA

MARCH™ Color-coded trauma treatment system		Mojo Crisis Responder Advanced	Mojo Crisis Responder Basic
	Product Number:	87-565-RC/BC	87-561-RC/BC
	Components	QTY	QTY
	Celox Rapid	2	
	Celox B-Con		2
	TMT Tourniquet	2	2
	Pre-Lubricated Nasal Airway	1	
	Positional Airway Card		1
	Sentinel Chest Seal	2	
	FoxSeal		1
	Compression Dressing		2
	Battle Bandage	2	
	Battle Wrap	2	2
	Rolled Gauze	2	2
	Emergency Blanket	2	2
	Mojo Crisis Responder Bag	1	1
	Pair Exam Gloves	2	2
	Trauma Shears	1	1
	MARCH QuickStart Guide	1	1