



## **MOJO® DIRECT ACTION AID BAG**

**From: \$507.06**

The Best Bring Their mojo

---

**SKU: 84-495-01-CB**

### **PRODUCT DESCRIPTION**

The Mojo® Direct Action Aid Bag strikes the balance between smaller Combat Lifesaver bags and bulkier advanced care bags. This ultra low profile, multi-mission medical bag is designed for use in close quarters and vehicle environments. The pockets, flaps, pouches and multi-loop retention system allow for maximum configuration options for a variety of patrol, assault, rescue and vehicle-based missions.

Inside the Direct Action Aid Bag are clear, removable pockets that allow for quick, organized treatments. And, if your mission should extend, the Direct Action Aid Bag is so versatile it can be adjusted quickly by adding the removable padded waist belt and three large external pockets to increase the amount of medical or mission gear stored. The Direct Action Aid Bag offers the field Medic immediate options for tailoring the kit to meet their ever changing operational needs.

Born from the real-life experiences of military special operations medics across years of deployments, “Mojo® Medic” Bags aim to simplify tactical medicine, even for the most advanced practitioners. Every pouch, pocket, strap and handle has been designed and purposefully placed to give first medics the best access to the life-saving equipment inside.

### **PRODUCT INFO**

- Ultra Low Profile
- Maximum Configurability
- Versatile Portability

## **Product Attributes**








- Mission configurable: organize your medical bag according to operational need
- Low profile allows for easy vehicle travel
- Multi-Loop retention system maximizes internal configurability
- Four large, clear front movable and/or removable procedure pockets
- Generous Velcro closure flaps on procedure pockets
- Removable padded waist belt

## **Clinical Benefits**

- External access hydration pocket
- Three large external removable equipment pockets for surgical operational equipment or extra medical supplies

## **SPECIFICATIONS**

- Unit Weight (stocked): 11.9 lbs
- Unit Dimensions: 24"L x 17"W x 7"D
- Berry Amendment & trade compliant
- Sourced, manufactured and assembled in the U.S.A.

MARCH™ Color-coded trauma treatment system		Mojo Direct Action Advanced	Mojo Direct Action Intermediate	Mojo Direct Action Basic
	Product Number:	84-497	84-493	84-495-01
	Components	QTY	QTY	QTY
	Celox Rapid	4	4	4
	TMT Tourniquet	4	4	4
	ITClamp	2		
	Prelubricated Nasal Airway	2	2	2
	Contol-Cric	1	1	
	Sentinel Chest Seal	4	4	2
	Dart			2
	Chest Decompression Set	2	2	
	Battle Bandage	3	3	
	Battle Wrap	2	2	
	Rolled Gauze	3	3	3
	Abdominal Dressing	1	1	1
	Compression Dressing			3
	Elastic Bandage 6"	2	2	2
	Elastic Bandage 4"			1
	Mojo Saline Lock Kit	1	1	1
	T.A.L.O.N. 7 Site IO Starter Set	1	1	
	Polycarbonate Eye Shield w/ Garter	2		
	Aluminum Eye Shield w/ Garter		2	2
	APLS LifeGuard Military	1		
	HAWK Warming Grid	1	1	
	Blizzard Light		1	1
	Mojo Direct Action Bag	1	1	1
	Mojo Cravat	2	2	2
	Universal Splint	2	2	2
	Pair Exam Gloves	8	8	8
	Trauma Shears	1	1	1
	Shears Retractor	1	1	1
	TCCC Card DD 1380	4	4	4
	Permanent Marker	2	2	2
	Cloth Medical Tape 3" x 10 yd	1	1	1
	Sharps Shuttle	1	1	1