



## **MOJO® HOLSTER BLEEDER KIT**

**\$89.23**

Stop Bleeding Even Faster

---

**SKU:** 88-115-03

### **PRODUCT DESCRIPTION**

The Mojo® Holster Bleeder Kit is an advanced, compact kit in a nylon pouch designed for the treatment of massive hemorrhage; the leading cause of preventable death following traumatic injury. This kit holds the TMT™ Tourniquet and Celox™ Rapid, which are easily deployed with the attached red, nylon pull strap. The pouch has a vertical belt loop and attaches directly to MOLLE. It can also be worn on a standard 2" wide utility and/or gun belt. The Mojo® Holster Bleeder Kit is available in both black and olive drab.

### **PRODUCT INFO**

- Treats Massive Hemorrhage, one of the five leading causes of traumatic death
- Ergonomically Designed for Ease of Use
- 5 Year Shelf Life

### **SPECIFICATIONS**

- Specifications:
  - Weight: 6.3 oz
  - Dimensions: 3.25" x 6" x 2.75"
- Single-Use
- 5 Year Shelf Life
- Berry Compliant, 500D nylon pouch
- 100% TAA Compliant