



MOJO® IFAK

From: \$50.25

Stand out with MARCH™ Inside

SKU: 81-132-01

PRODUCT DESCRIPTION

Mojo® IFAK uses MARCH™ medical gear in a user-friendly, color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five preventable causes of death in combat.

The Mojo® IFAK is configurable and extremely versatile for every mission possibility. It can be worn in a drop leg configuration, and with the removable leg strap it can also attach directly to MOLLE vest/harness, or to a 2" duty/gun belt. The Mojo® IFAK is made of 1000D Cordura® and available with three different content configurations, as well as four color options.

- Bag Only- PN 81-130
- Basic- PN 81-132-01
- Intermediate- PN 81-133
- Advanced- PN 81-135

Image shown is PN 81-135, actual product may vary

PRODUCT INFO

- Quick Access
- Intuitive Layout
- Versatile Portability

Product Attributes

- Ergonomic design evenly distributes weight and is worn comfortably on your body during long missions








- Wearable in two ways (belt, MOLLE system)
- Made of 1000D Cordura®
- Systematic access and full visualization of your life saving equipment ensures the smoothest casualty intervention
- 180-degree top-opening zipper panel for quick grab access
- Fully mission-configurable kit

Clinical Benefits

Mojo® IFAK Basic uses MARCH™ medical gear in a user-friendly, color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five preventable causes of death in combat.

SPECIFICATIONS

- PN: 81-135
 - coyote brown, black, foliage green, multicam
- Up to 5-Year Shelf Life
- Single Use
- Specifications:
 - Dimensions: 6.2" x 5.7" x 3.3"
 - Weight: 1.75 lbs
- Berry Amendment / Trade Compliant
- Made In USA
- Made of 1000D Cordura

MARCH™ Color-coded trauma treatment system		Mojo IFAK Advanced	Mojo IFAK Intermediate	Mojo IFAK Basic
	Product Number:	81-135	81-133	81-132-01
	Components	QTY	QTY	QTY
	Celox Rapid	1	1	1
	TMT Tourniquet	1	1	1
	ITClamp	1		
	Prelubricated Nasal Airway	1	1	1
	Sentinel Chest Seal	1	1	1
	Chest Decompression Set	1	1	
	Compression Dressing			1
	Battle Bandage	1	1	
	Battle Wrap	1	1	
	Rolled Gauze	1	1	1
	Polycarbonate Eye Shield w/ Garter	1		
	Aluminum Eye Shield w/ Garter		1	1
	Mojo IFAK Bag	1	1	1
	Pair Exam Gloves	1	1	1
	Trauma Shears	1	1	
	Shears Retractor	1	1	
	TCCC Card DD 1380	1	1	1
	Permanent Marker	1	1	1
	MARCH QuickStart Card Set	1	1	1