

**MOJO® PATROL IFAK****From: \$71.50**

Stand out with MARCH™ Inside

---

**SKU: 81-121****PRODUCT DESCRIPTION**

The Mojo® Patrol IFAK is specifically designed for lightweight and low-profile wear. It has a quick-release handle for single-handed deployment that can be stowed when not needed. The Mojo® Patrol IFAK fully opens for extremely generous access and can attach to a 2" duty belt or directly to MOLLE. It uses MARCH™ medical gear in a color-coded trauma treatment system to treat the five leading preventable causes of death on the battlefield. The Patrol IFAK is made of 1000D Cordura® nylon, and is available in four different color options, as well as three different content configurations.

- Bag Only- PN 81-120
- Basic- PN 81-121
- Intermediate- PN 81-123
- Advanced- PN 81-125

\*Image shown is PN 81-125, actual product may vary\*

## PRODUCT INFO

- Quick Access
- Intuitive Layout
- Versatile Portability

## Product Attributes








- Ergonomic design evenly distributes weight and is worn comfortably on your body during long missions
- Wearable in two ways (belt, MOLLE system)
- Made of 1000D Cordura®
- Systematic access and full visualization of your life saving equipment ensures the smoothest casualty intervention
- 180-degree top-opening zipper panel for quick grab access
- Fully mission-configurable kit

## Clinical Benefits

Mojo® Patrol IFAK Basic uses MARCH™ medical gear in a user-friendly, color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five preventable causes of death in combat.

## SPECIFICATIONS

- NSN: Pending
- PN: 81-125
  - (coyote brown, black, foliage green, multicam)
- Up to 5 year shelf life
- Single Use
- Specifications:
  - Weight: 1.7 lbs
  - Dimensions: 7.8 x 6.1 x 3.3
- Available Direct, Prime Vendor, ECAT and GSA
- Berry Compliant / Trade Compliant
- Made in USA

MARCH™ Color-coded trauma treatment system		Mojo Patrol IFAK Advanced	Mojo Patrol IFAK Intermediate	Mojo Patrol IFAK Basic
	Product Number:	81-125	81-123	81-121
	Components	QTY	QTY	QTY
	Celox Rapid	1	1	
	TMT Tourniquet	1	1	1
	Prelubricated Nasal Airway	1	1	
	Sentinel Chest Seal	1		
	Dart	1		
	FoxSeal		1	1
	Compression Dressing		1	1
	Battle Bandage	1		
	Battle Wrap	1		
	Rolled Gauze	1	1	1
	Polycarbonate Eye Shield w/ Garter	1		
	Aluminum Eye Shield w/ Garter		1	1
	Mojo Patrol IFAK Bag	1	1	1
	Pair Exam Gloves	1	1	1
	TCCC Card DD 1380	1	1	
	Permanent Marker	1	1	
	MARCH QuickStart Card Set	1	1	1