



MOJO® POCKET BLEEDER KIT

From: \$60.45

Massive Hemorrhage Controlled

SKU: 88-100

PRODUCT DESCRIPTION

The Mojo® Pocket Bleeder Kits are designed to treat massive hemorrhage; the leading cause of preventable death from traumatic injury. These kits come in a compact, vacuum sealed, nylon pouch, which can be easily carried in any pocket. The Pocket Bleeder Kits use MARCH™ medical gear, which is a color-coded, trauma treatment system distinctly packaged with icon identifiers. All Mojo® Pocket Bleeder Kits include a tourniquet, Celox™ Rapid hemostatic gauze, pair of exam gloves, and MARCH™ Quick Reference Card Set. These convenient kits are offered in four different content configurations to accommodate any departmental needs.

PRODUCT INFO

- Treats Massive Hemorrhage
- Compact and Easy to Use
- Superior Performance

SPECIFICATIONS

- Weight: 5.2 oz -6.9 oz
- Dimensions: 3.5" x 6" x (1.5"-3.5", widths vary based on configuration)
- Single-Use
- Not made with natural rubber latex
- 5 year shelf life
- 100% TAA Compliant
- Berry compliant nylon carriers

