



## **MOJO® PUBLIC ACCESS TRAUMA STATION**

**From: \$1,065.90**

Empowering the Public to Save Lives

---

**SKU: 87-600-01**

### **PRODUCT DESCRIPTION**

The Mojo® Public Access Trauma Station contains everything needed to treat the five leading causes of preventable death in traumatic injuries. The kit contains a nylon carrier inside a wall mounted box, making life saving interventions readily available. It is recommended this kit is placed in common areas similar to public access Automated External Defibrillators.

Using only MARCH™ medical gear allows for simplified training that visually cues the user to the right product, the right intervention at the right time. Includes Celox® RAPID hemostatic gauze, which is clinically effective in the presence of blood thinners<sup>1</sup>; and the versatile TMT tourniquet which accommodates varying limb circumferences from Minimum 5 3/4" to Maximum 35 1/4".

### **PRODUCT INFO**

#### **Capabilities**

- Treats 6-12 Victims
- Wall Mounted Box Included
- Quick Start Guide

#### **MARCH Public Access Trauma Station (Advanced, Intermediate, Basic) Contain:**

TMT Tourniquet (2)

CPR Microshield (2)

Exam Gloves (4 pr)

Trauma Shears (2)

Positional Airway Instruction Card (2)

MARCH Quick Reference Cards (2)

Public Access Nylon Carrier (1)

Wall Mount Case (1)

Individual Trauma Kit (6) - (available in Advanced, Intermediate, or Basic):

**In addition to the above, the stations also contain the below:**

Advanced Individual Trauma kit contains:

Celox Rapid (1)

TMT Tourniquet (1)

Sentinel Chest Seal (2)

Emergency Blanket (1)

Exam Gloves (1 pr)

Rolled Gauze (1)

4 "Compression Dressing (1)

Positional Airway Instruction Card (1)

MARCH Quick Reference Guide (1)

Intermediate Individual Trauma Kit contains:

Celox Rapid (1)

SWAT-Tourniquet (1)

Sentinel Chest Seal (2)

Emergency Blanket (1)

Exam Gloves (1 pr)

Rolled Gauze (1)

4 "Compression Dressing (1)

Positional Airway Instruction Card (1)

MARCH Quick Reference Guide (1)

Basic Individual Trauma Kit contains:

Celox Rapid (1)

SWAT-Tourniquet (1)

FoxSeal (1)

Emergency Blanket (1)

Exam Gloves (1 pr)

Rolled Gauze (1)

4 "Compression Dressing (1)

Positional Airway Instruction Card (1)

MARCH Quick Reference Guide (1)

## **Clinical Benefits**

- Each station addresses the five leading causes of death following trauma
- Contains tourniquets that accommodate varying limb circumferences
  - Min: 5 3/4"
  - Max: 35 1/4"
- Hemostatic gauze is proven to work in therapeutically anti-coagulated patients
- Color-coded packaging for easy identification during stressful situations

## **SPECIFICATIONS**

- Nylon Carrier Color- Red
- Case Dimensions:
  - 15"L x 11.5" W x 6" D
- Three configurations:
  - 87-600-01 (Basic)
  - 87-600-02 (Intermediate)
  - 87-600-03 (Advanced)
- TAA Compliant
- Assembled in USA

## REFERENCES

- 1 Hemostatic effect of a chitosan linear polymer (Celox®) in a severe femoral artery bleeding rat model under hypothermic or warfarin therapy. Koksai O et al Turk J Trauma & Emerg Surg. 2011; 17 (3):199-204