








MOJO® RESPONDER BELTPACK**From: \$100.10**

SKU: N/A**PRODUCT DESCRIPTION**

The Responder Beltpack is a tactically stocked trauma kit. Its mid-size and elongated belt design, straps medical gear around the waist so they can be easily accessible during stressful situations. The Responder Beltpack is intended for first responders and rescue task force members in emergency situations. It uses MARCH™ medical gear with icon-identifiers, in a color-coded treatment system targeting the five leading causes of preventable death following trauma. The Responder Beltpack is offered in two different configurations to accommodate the needs of any department or situation.

SPECIFICATIONS

- Weight: 3 lbs
- Dimensions: 15" L X 7" W X 4.75" H
- Available in red or black diamond, check pattern
- Up to 5 year shelf life
- Made in the USA

| MARCH™ Color-coded trauma treatment system | | Responder Beltpack Advanced | Responder Beltpack Basic |
|---|-----------------------------|--------------------------------|-----------------------------|
|  | Product Number: | 87-575-RC/BC | 87-571-RC/BC |
| | Components | QTY | QTY |
|  | Celox Rapid | 3 | |
| | Celox B-Con | | 3 |
| | TMT Tourniquet | 3 | 1 |
| | SWAT-T | | 2 |
|  | Pre-Lubricated Nasal Airway | 2 | |
| | Positional Airway Card | | 1 |
|  | Sentinel Chest Seal | 2 | |
| | FoxSeal | | 1 |
|  | Compression Dressing | | 2 |
| | Battle Bandage | 3 | |
| | Battle Wrap | 3 | 3 |
| | Rolled Gauze | 2 | 3 |
|  | Emergency Blanket | 2 | 2 |
|  | Responder Beltpack Carrier | 1 | 1 |
| | Pair Exam Gloves | 4 | 4 |
| | Trauma Shears | 1 | 1 |
| | Shears Retractor | 1 | |
| | MARCH QuickStart Guide | 1 | 1 |