



MOJO® SALINE LOCK KIT

\$12.99

SKU: 80-110

PRODUCT DESCRIPTION

The Mojo® Saline Lock Kit allows for quick and effective extremity venous cannulation under the most demanding conditions in one complete set. It contains the necessary items for starting and securing a saline lock. The rugged, low-cube, vacuum-sealed pouch keeps the contents of the kit in prime condition until needed.

PRODUCT INFO

- PN: 80-110
- Shelf Life: 6 - 18 months
- Single Patient Use
- Not made with natural rubber latex
- Specifications:
 - Dimensions: 3"W x 6.5"H x 1"D
 - Weight: 1.7 oz
- Available: ECAT, PV, VIPA, GSA, CEC
- Made in USA

SPECIFICATIONS

Mojo® Sakine Lock Kits

End Item:	Basic
Product Number:	88-110
Components:	QTY
18 GA Catheter	1
20 GA Catheter	1
Saline Flush (5cc pre-filled syringe)	1
Povidone Pad	1
Alcohol Prep Pad	2
Needle-Free Valve w/ Luer Lock	1
Constricting Band	1
Transparent Film Dressing	1
2"x2" Gauze Pad (2 pack)	1