



MOJO® VEHICLE RESCUE

From: \$577.00

Rescue that goes where you go

SKU: 85-200








Info:: Unfortunately, due to supply chain challenges, all kits that contain the iTClamp will ship without it. This product will be shipped to you once available at no additional cost. Please contact us directly at 855-428-6074 opt #2 with any questions. |

PRODUCT DESCRIPTION

The Mojo® Vehicle Rescue is a rapidly deployable kit that keeps evacuation easy. This Mojo system includes the Mantis™ tri-fold evacuation board, that is standardized in the USMC and SOCOM kits with the Mojo® CLS Intermediate Bag attached to the front. It is ideal for quick vehicle extraction due to its simple, rapid deployment. The Mojo® Vehicle Rescue uses MARCH™ medical gear in a user-friendly, color-coded trauma treatment system distinctly packaged with icon identifiers to treat the five leading causes of preventable death. Mojo® Vehicle Rescue (Intermediate) comes complete with the Mantis Evacuation Board, CSI Bag®, AllEvac® Patient Restraint System, and Mojo® CLS Intermediate Kit. It is made of 1000D Cordura, and is available in four different colors options.

(Intermediate Bag Pictured)

PRODUCT INFO

MARCH™ Color-coded trauma treatment system		Mojo Vehicle Rescue Advanced	Mojo Vehicle Rescue Intermediate	Mojo Vehicle Rescue Basic
	Product Number:	85-205	85-203	85-201
	Components	QTY	QTY	QTY
	Celox Rapid	2	2	2
	TMT Tourniquet	2	2	2
	ITClamp	2		
	Prelubricated Nasal Airway	2	1	1
	Sentinel Chest Seal	2	2	2
	Dart			2
	Chest Decompression Set	2	2	
	Compression Dressing			2
	Elastic Bandage 6" x 4.5 yd			1
	Abdominal Dressing	1	1	1
	Battle Bandage	2	2	
	Battle Wrap	2	1	
	Rolled Gauze	2	2	2
	Polycarbonate Eye Shield w/ Garter	2		
	Aluminum Eye Shield w/ Garter		1	1
	Blizzard Light Blanket		1	1
	HAWK Warming Grid	1	1	
	APLS LifeGuard Military	1		
	Mojo CLS Carrier	1	1	1
	AllEvac Mantis System	1	1	1
	Mojo Cravat	3	3	3
	Hook Knife	1	1	
	Pair Exam Gloves	4	4	4
	Universal Splint	1	1	1
	Cloth Medical Tape 3" x 10 yd	1	1	1
	Trauma Shears	1	1	1
	Shears Retractor	2	2	
	TCCC Card DD 1380	2	2	2
	Permanent Marker	2	2	2
	MARCH QuickStart Guide	1	1	1

SPECIFICATIONS

- PN: 85-203 (four color options available)
 - coyote brown, black, foliage green, multicam
- Up to 5-Year Shelf Life
- 500 lb load rating
- Specifications (Kitted):
 - Dimensions: 26.5" x 17.5" x 5.5" D
 - Weight: 25 lbs
- Specifications (Mantis Evacuation Board):

- Folded Dimensions: 25" L x 16" W x 3.5" D
- Deployed Dimensions: 72" L x 16" W x 1" D
- Weight: 12.8 lbs
- Berry Compliant
- Made In USA