



NIO ADULT

\$125.00

The NIO Adult is an automatic intraosseous access device packaged for quick, safe, easy vascular access in **10 seconds or less***. No batteries. No assembly. No extra parts.

SKU: NIO-A

NSN: 6515-01-667-2884 |

PRODUCT DESCRIPTION

When time is life, trust the NIO Adult™ to deliver safe, reliable vascular access in 10 seconds¹ — no batteries, assembly, or extra parts required. Designed with providers in mind, the preassembled NIO Adult™ features no exposed needle, virtually eliminating the risk of needlestick. With first attempt success rates as high as 98%², the NIO Adult™ supports common user errors, empowering providers to deliver safe and effective patient care when it matters most. The NIO Adult™ accommodates adult patients 12 years of age and up.

1. Deployment of the NIO - Next-Generation IO™ Device by Medical Personnel in a Cadaver Lab. 2020.
2. Swisa, Einat. NIO-A.V. Clinical Studies Report, NIO Human Cadaver Studies. 2013.

PRODUCT INFO

- Vascular access in 10 seconds or less¹
- Automatic device, no batteries, assembly, or extra parts required
- Pocket-sized and lightweight
- Sterile and single-use device—eliminates the risk of cross-contamination
- Self-contained device with no exposed needle
- Pre-set needle penetration depth
- 5-year shelf life (at room temp.)
- Hyper-realistic training system
- Eliminates the risk of CLABSI

1. Deployment of the NIO - Next-Generation IO™ Device by Medical Personnel in a Cadaver Lab.

2020.

SPECIFICATIONS

- PN: NIO-A
- NSN: 6515-01-667-2884
- Needle gauge: 15G
- Device dimensions: 5.47in x 2.24in x 1.07in (13.9cm x 5.7cm x 2.7cm)
- Package dimensions: 7.64in x 3.3in x 1.3in (19.4cm x 8.4cm x 3.3cm)
- Weight: 3.2oz (91g)
- Shelf life: 5 years at room temperature