



NIO+ ADULT

From: \$189.00

With the NIO+ Adult, providers can achieve safe, reliable vascular access in 10 seconds

SKU: NIO+

PRODUCT DESCRIPTION

The NIO+ Adult™ is a durable, reliable IO device built to withstand harsh conditions and challenging environments. Packaged in an ultra-compact hard-shell case that easily stores in a small medical pouch, the NIO+ Adult™ is a practical solution for limited-storage scenarios. The NIO+ Adult™ features a reinforced core, which resists elevated thermal temperatures, making it an ideal solution for extreme environments.

PRODUCT INFO

- Tough, rapid, and reliable: the NIO+ Adult™ withstands harsh climates and challenging environments to deliver IO access in 10 seconds¹.
- 5-year extreme temperature shelf life*: The NIO+ Adult™ can be stored either at room temperature for 5 years or at 122°F (50°C) for up to 5 years in the US (up to 3 years in EU countries), making it an ideal solution for extreme environments.
- Ultra-compact hard-shell case that easily stores in a small medical pouch
- Automatic device, no batteries, assembly, or extra parts required
- Sterile and single-use device—eliminates the risk of cross-contamination
- Self-contained device with no exposed needle
- Pre-set needle penetration depth

* Deployment of the NIO - Next-Generation IO™ Device by Medical Personnel in a Cadaver Lab. 2020.

SPECIFICATIONS

- **PN:** NIO+
- **Needle gauge:** 15G
- **Device dimensions:** 5.47in x 2.24in x 1.07in (13.9cm x 5.7cm x 2.7cm)
- **Package dimensions:** 6.5in x 2.8in x 1.8in (16.6cm x 7.2cm x 4.6cm)

- **Weight:** 5.1oz (145.5g)
- **Shelf life:** 5 years at room temperature
 - US: up to 5 years at 122°F (50°C)
 - EU: up to 3 years at 122°F (50°C)