



## **POCKET MASK**

**\$8.00**

Protecting the Lifesaver

---

**SKU:** 23-025

**NSN:** 6515-01-215-4177 |

## **PRODUCT DESCRIPTION**

The Pocket Mask is designed for mouth to mask ventilation of a non-breathing adult. It is also used as a barrier that will direct expired air from the patient away from the user. The Pocket Mask is equipped with a filter-protected, replaceable, one-way valve. It folds flat for a more compact size and is stored in a splash-resistant, clam shell storage case. Alcohol wipe and one pair of nitrile gloves are included. This device is intended to be used by persons who have received adequate training.

## **PRODUCT INFO**

- High efficiency and low resistance one-way valve minimizes the possibility of cross-contamination
- 15mm O.D. connector can be attached to any standard ventilation and oxygen equipment
- Ergonomically shaped mask with soft-air cushion for leak-free performance, easy grasp, and patient comfort
- Disposable filter to avoid cross-contamination
- Elastic green strap to keep mask in place on patient
- Mask is reusable and easy to clean
- Instructions for use, disposable gloves, alcohol wipe, and pocket storage case also included

## **SPECIFICATIONS**

- NSN: 6515-01-215-4177
- PN: 23-025
- Single Patient Use
- Not made with natural rubber latex

- Operational Temperature: -18° C to 40° C
- Specifications:
  - Dimensions: 5.1" x 4.2" x 1.9"
  - Weight: 4.5 oz
- Available Direct, Prime Vendor, ECAT
- 100% TAA Compliant