



## **POLYCARB EYE SHIELD**

**\$5.50**

Strong and Lightweight Eye Protection

---

**SKU:** 06-205

## **PRODUCT DESCRIPTION**

The MARCH™ Polycarb Eye Shield is a standard ophthalmic dressing that meets current military pre-hospital guidelines for the treatment of traumatic eye injury in the field. It is a polycarbonate eye shield that is enhanced with an attached edge protecting garter. These eye shield units not only save time, but they are more comfortable and better tolerated by casualties.

## **PRODUCT INFO**

### **Capabilities**

- Protective
- Comfortable
- Reliable

### **Product Attributes**

- Polycarbonate formed eye shield
- Comfortable fabric garter around edge
- Strong and lightweight

### **Clinical Benefits**

- Blocks external pressure to injury site
- Accommodates either eye
- Meets CoTCCC treatment protocols for penetrating eye trauma

## **SPECIFICATIONS**

- NSN: None
- PN: 06-205
- Specifications
  - Unit Dimensions
    - Packaged 2.5" x 3.25" x 0.4"
  - Unit Weight: 0.3 oz
- Available Direct, Prime Vendor, ECAT and GSA
- 100% TAA Compliant
- Made in USA