

PREVAHEXCHX™

\$234.25

Quantity: Box of 50

SKU: 80-812

Info:: Quantity: Box of 50 |



PRODUCT DESCRIPTION

PrevaHexCHX™ is a thin, transparent, antimicrobial film dressing. The dressing incorporates groundbreaking advanced technology utilizing pure chlorhexidine, providing rapid, broad-spectrum protection of wound and catheter sites. PrevaHexCHX is chemically composed of a polyurethane film with an acrylic pressure sensitive adhesive film. Protects IV sites and is a component of our Valkyrie LTOWB Sets.

SPECIFICATIONS

- Part Number: 80-812
- Dimensions (EA): 4”L x 4.75” W
- Weight (EA): 0.3 oz
- Sterile, Single-use
- Rapid Onset of Action
- Free of PCA (p-chloroaniline)
- Effective Against Multiple Drug Resistant Microbes
- No Skin Irritating Salts, Acids, Binders, or Colophony
- Shelf Life: 3 Years
- Trade Compliant
- Made in the USA