

**PRO PELVIC SPLINT™****\$56.00**

Pelvic Immobilization

---

**SKU: 55-110****PRODUCT DESCRIPTION**

The Pro Pelvic Splint™ is designed for emergency management of suspected pelvic fractures. It is extremely lightweight, simple to apply, and can be adjusted by a single operator. Using the preferred method of application, the neoprene band can be trimmed once the splint has been sized to the patient. This allows greater access to the inguinal region and urethra for urinary catheterization and/or femoral vein/artery access. The Pro Pelvic Splint is suitable for all size patients.

**PRODUCT INFO**

- Made of 6mm Neoprene fabric
- One size fits all
- May be cut to size to allow access to abdomen and groin without compromise
- Soft material with no stitching reduces risk of pressure ulceration and/or friction
- The perineum is left uncovered

**SPECIFICATIONS**

- PN: 55-110
- Specifications:
  - Packaged Dimensions: 12" x 8" x 1.5"
  - Deployed Dimensions: 7.5" x 48" x 0.25"
  - Unit Weight: 9.5 oz
- Single Patient Use
- Latex Free
- 100% TAA Compliant

