

## ROLLED GAUZE

**\$3.77**

What are You Packing?

---

**SKU:** 33-300



## PRODUCT DESCRIPTION

The MARCH™ Rolled Gauze is a sterile and compact 6-ply; 4.5" x 4.1yd gauze used for the treatment of bleeding and bandaging wounds. This vacuumed-sealed gauze is made of high quality U.S. cotton and packaged in a versatile configuration that is ideal for both linear wound packing and basic bandaging and wrapping. Gauze is an essential item for any medical bag and can be used to control hemorrhage in conjunction with a compression bandage, used as backing gauze for hemostatic dressings or for wound bandaging.

## PRODUCT INFO

### Capabilities

- Configured for wound packing and wrapping
- Easy-To-Open tear notches
- Ideal Adjunct to compression bandage

### Product Attributes

- 4.5" x 4.1 yd Sterile, 6-ply cotton gauze
- Made in the U.S. & Berry Amendment Compliant

### Clinical Benefits

- Perfect for wound bandaging and wrapping when using exterior of roll
- Linear wound packing when pulling from interior of roll
- Cotton provides maximum absorption to help control bleeding

- Expands to cover larger wound areas
- Ideal backing for hemostatic agents

## **SPECIFICATIONS**

- PN: 33-300
- Sterile / Single-Use
- 6-Year Shelf Life
- Specifications:
  - Unit Dimensions
    - Packaged: 6" x 2" x 2"
    - Deployed: 4.5" x 4.1yd
  - Unit Weight: 1.5 oz
- Available Direct, Prime Vendor, ECAT and GSA
- 100% TAA Compliant
- Made in USA

## **TRAINING**

[Rolled Gauze Label & IFU](#)