

**SPIDER STRAP XL®****\$183.75**

Ultimate Patient Immobilization

**SKU: N/A****PRODUCT DESCRIPTION**

The Spider Strap XL is a patient immobilization system. It is made of polypropylene with visible, light reflective strips. It also has webbing, with hook and loop fasteners, that help to secure the patient. The Spider Strap XL has 12 points of attachment compared to the traditional 10. The cross straps are longer to accommodate any gear. A storage carrier completes the set.

**PRODUCT INFO**

- Includes carrying pouch
- No assembly required
- Lightweight and durable
- Visible, light-reflective tape on attachment points
- Secures and stabilizes
- Simple instructions for use
- 12 Points of attachment for secure placement
- Indefinite shelf life

**SPECIFICATIONS**

- PN: 53-589-OD (Olive Drab)
  - NSN: 6530-01-593-0010
- PN: 53-589-BL (Black)
  - NSN: 6530-01-586-2411
- Specifications:
  - Packaged Dimensions: 18.5" x 4"

- Unit Weight: 14 oz
- Cycle PSI Average: 200 lbs.
- Shear Strength at Length: 220 lbs.
- Shear Strength Width: 170 lbs.
- Berry Compliant
- Made in USA

**NSNS****PN:** 53-589-BL**NSN:** 6530-01-586-2411**PN:** 53-589-OD**NSN:** 6530-01-593-0010