

**T.A.L.O.N.™****\$2,022.25**

7-Site Tactical Advantage

SKU: 35-081**NSN:** 6515-01-642-7919 | **Info::** Quantity:
Box of 6 |

PRODUCT DESCRIPTION

The T.A.L.O.N.™ (Tactically Advanced Lifesaving IO Needle) is the single MARCH™ system solution for intraosseous (IO) vascular access. The needle allows providers rapid and accurate IO access to all seven sites. The TALON delivers life-saving fluids and medications rapidly and accurately for casualty resuscitation.

T.A.L.O.N.™ is the first and only FDA-cleared manual IO device with 7-site access: sternum, bilateral proximal humerus, bilateral proximal tibia, and bilateral distal tibia. A single needle set provides IO access—there are no powered tools or multiple needles necessary and no special removal tool is required. Simple step-by-step instructions and minimal insertion pressure is required to ensure that combat lifesavers and medics can use the needle set effectively with proper training.

PRODUCT INFO

- **7-Site IO Vascular Access**
- **Simple Delivery** of life-saving fluids and medication
- **Rapid Resuscitation**
- Access to any IO site with one manual device
- No need for additional devices, multiple needles or special removal tools
- Meets CoTCCC recommendations for IO devices
- First and only FDA-cleared manual IO device with 7-site access:
 - Sternum
 - Bilateral proximal humerus

- Bilateral proximal tibia
- Bilateral distal tibia
- Delivers fluids and medications rapidly and safely

SPECIFICATIONS

- NSN: 6515-01-642-7919
- PN: 35-081
- FDA 510(k) cleared
- Sterile / Single-Use
- 4-year shelf life
- Specifications:
 - Unit Dimensions
 - Packaged: 8.25"x 5.1"x 2"
 - Unit Weight: 1.4oz
 - Not made with natural rubber latex
- Available Direct, Prime Vendor, ECAT, CEC and GSA
- 100% TAA Compliant
- Made in USA

TRAINING

[T.A.L.O.N. Video](#)