



TMT™ RIGID HOLSTER WITH DCL ATTACHMENT

From: \$37.00

SKU: N/A

PRODUCT DESCRIPTION

The TMT Rigid Holster with DCL Combat Loop Belt Attachment is a custom designed Kydex® holster for the CoTCCC Recommended Tactical Mechanical Tourniquet (TMT™). The DCL mounting attachment enables the holster to be carried in the vertical or horizontal position on belts and web gear. This attachment works like a door swinging closed and snapping into place for a secured carry. The DCL can be unlocked and removed via a one handed spring closure release. This same release has a built in slide lock safety to secure the mounting to your belt. The DCL hole spacing features 3/4" inch centers and has a removable belt size adjustment bar that fits 2" inch belts down to 0.25" inch cords and webbing.

PRODUCT INFO

- Holster constructed of Kydex ®
- Adjustable sizing
- Adjustable friction screws
- 360 swivel
- Waterproof, Washable
- Scratch resistant
- Rugged, Secure
- Lightweight
- Holds its shape better; will not stretch or shrink in normal conditions

SPECIFICATIONS

- PN: 89-121-XX (Empty Holster)
- PN 89-321-XX (TMT Holster Set)
- Latex Free
- Non Shelf Life Product
- Specifications:
 - Weight:
 - Empty Holster: 3.2 oz
 - TMT Holster Set: 6.1 oz
 - Dimensions: 3.4"L x 3.5" W x 2"D
- Assembled in USA
- Available in Black, Coyote Brown, Olive Drab, Multicam (*Multicam pricing may vary)