



TMT™ RIGID HOLSTER WITH MALICE CLIP® ATTACHMENTS

From: \$37.00

SKU: N/A

PRODUCT DESCRIPTION

The TMT Rigid Holster with Malice Clip Attachments is a custom designed Kydex® holster for the CoTCCC recommended Tactical Mechanical Tourniquet (TMT™). Malice Clips are injection molded, heavy duty, reusable connecting clips that are designed to attach nylon carriers to MOLLE/PALS style equipment. Malice Clips are not affected by heat or cold, and will never corrode or lose their subdued finish. Once properly closed, the clip will not open until it is disengaged by the user using a flat tipped object such as a knife or screwdriver.

PRODUCT INFO

- Holster constructed of Kydex ®
- Injection molded connection clips
- Malice Clip attaches to MOLLE
- Adjustable friction screws
- Waterproof, Washable
- Scratch resistant
- Rugged, Secure
- Lightweight
- Holds its shape better; will not stretch or shrink in normal conditions

SPECIFICATIONS

- PN: 89-123-XX (Empty Holster)
- PN 89-323-XX (TMT Holster Set)
- Latex Free

- Non Shelf Life Product
- Specifications:
 - Weight: 3.2 oz
 - Dimensions: 3.4"L x 3.5" W x 2"D
- Made in USA
- Available in Black, Coyote Brown, Olive Drab, Multicam (*Multicam pricing may vary)